

S E T T O R E

laboratorio

laboratory
balances



0111111111

LABORATORY ELECTRONIC BALANCES & OENOLOGICAL EQUIPMENTS

ISO 9002



SIT n. 94



Carl Lidén Vågservice
Huskvarnavägen 40
554 54 Jönköping

Tel: 036-12 57 80
Fax: 036-12 57 89

sales@carlliden.com
<http://www.carlliden.com>

**SEMIMICRO ANALYTICAL BALANCE DUAL RANGE
mod. E/50S**

Semimicro:

- Capacity: 90 g
- Readability: ± 0.01 mg
- Linearity: ± 0.02 mg
- Stabilisation time: 8/10 sec

Macro:

- Capacity: 202 g
- Readability: ± 0.1 mg
- Linearity: ± 0.2 mg
- Stabilisation time: 5/8 sec

- Pan dimension: diam. 80 mm
- Built-in calibration weight
- Autocalibration
- Serial output RS232 I/O
- Facility for weighing underneath the balance
- Acid resistant membrane keyboards with MODE/PRINT, TARE-ZERO, ON/OFF keys
- Reading by 7-LEDs display
- Stable weight determination signal
- Bubble and feet for levelling
- Operating temperature: 10/35°C
- AC power requirements/voltage: 220V -15%/+10% (on request 110 V)
- Absorption: 15 VA
- Weighing chamber dimensions: W 190, D 160, H 250 mm
- Total dimension: W 200, D 410, H 300 mm
- Net weight 11 kg

On request:

legal metrologic approval (EEC Directive 90/384)

Carl Lidén Vågservice
Huskvarnavägen 40
554 54 Jönköping

Tel: 036-12 57 80
Fax: 036-12 57 89

sales@carlliden.com
<http://www.carlliden.com>