

S E T T O R E

laboratorio

laboratory
balances



IN VIVO DESIGN

LABORATORY ELECTRONIC BALANCES & OENOLOGICAL EQUIPMENTS

ISO 9002



SIT n. 94



Carl Lidén Vågservice
Huskvarnavägen 40
554 54 Jönköping

Tel: 036-12 57 80
Fax: 036-12 57 89

sales@carlliden.com
<http://www.carlliden.com>

Max 7500 g - d = 0,1 g

- STRENGTHENED ANTI-SHOCK SYSTEM
- PIECE COUNTING FUNCTION
- PRACTICALLY INDESTRUCTIBLE WEIGHING CELL WITH NO TRANSPORT PROBLEMS
- UNIVERSAL POWER SUPPLY AVAILABLE ALL OVER THE WORLD (110 - 240 V)
- IT STORES IN MEMORY THE WEIGHT IF THE ELECTRICITY SUPPLY IS LOST
- PAN 190 x 190 mm - LINEARITY $\pm 0,3$ g
- LIGHT (4,5 Kg) - LESS TRANSPORT CHARGES
- SERIAL OUTPUT RS232
- SUBJECT TO TECHNICAL CHANGES WITHOUT NOTICE.

Carl Lidén Vågservice
Huskvarnavägen 40
554 54 Jönköping

Tel: 036-12 57 80
Fax: 036-12 57 89

sales@carlliden.com
<http://www.carlliden.com>